Fermathron™

Manufacturer:

Hyaltech Ltd.
Heriot Watt Research Park
Edinburgh, EH14 4AP, UK
www.hyaltech.com

Distributor:

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Fermathron™

Sodium hyaluronate 20 mg/2.0 ml

For the relief of pain and stiffness of the knee joint and other synovial joints in patients with degenerative and traumatic changes to the synovial joint.

Presentation

Fermathron™ is a clear solution of sterile 1% sodium hyaluronate in a phosphate buffered saline contained in a pre-filled syringe for single intra-articular injection into the synovial space of the joint.

2.0 ml of Fermathron™, sterilised by filtration, is enclosed within a glass, ready to use, disposable syringe. The syringe is packed within a blister pack and an outer cardboard carton. The outer surface of the pre-filled syringe has been sterilised by ethylene oxide.

Dosage and Administration

The dosage regimen is injection into the affected synovial joint space once a week for up to five injections depending on the severity of the degenerative or traumatic change to the synovial joint.

The recommended dosage regimen for patients with mild to moderate osteoarthritis of the knee joint is up to five weekly injections of 2.0 ml into the synovial space of the knee joint.

Clean the skin around the injection site with antiseptic and allow to dry before the injection is given.

If joint effusion is present it should be aspirated before injection of Fermathron™.

The contents of the syringe are sterile and should be injected using a sterile needle of an appropriate size (19 to 20 gauge is recommended). The syringe is fitted with a Luer lock (6%).

Discard the syringe and needle after single use.

Uses

For the relief of pain and stiffness of the knee joint and other synovial joints in patients with degenerative and traumatic changes to the synovial joint. The duration of effect in patients with mild to moderate osteoarthritis is up to six months.

The performance of Fermathron™ is due to its biocompatibility and physicochemical properties. The sodium hyaluronate contained in Fermathron™ is a biopolymer composed of repeating disaccharide units of N-acetylglucosamine and glucuronic acid and though it is biosynthesised by the bacterium *Streptococcus equi* it has been shown to be the same as the sodium hyaluronate which is found in the human body. The hyaluronan supplements the hyaluronan found naturally in the synovium but which has been depleted by degenerative and traumatic changes to the synovial joint.

Contra-indications

Do not inject Fermathron™ if the area of the injection is infected or where there is evidence of skin disease.

Special warnings and precautions

Do not use if packaging has been damaged. Do not use after the expiry date.

Sodium hyaluronate is manufactured by fermentation of *Streptococcus equi* and rigorously purified. However, the physician should consider the immunological and other potential risks that can be associated with the injection of any biological material. Do not use for children.

Adverse reactions

Transient pain and swelling may occur with intra-articular injections.

Rarely, an inflammatory reaction could occur which may or may not be associated with Fermathron™.

Use during Pregnancy and Lactation

There is no evidence concerning the safety of Fermathron™ in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the orthopaedic surgeon.

Incompatibilities

Fermathron™ has not been tested for compatibility with other substances for intra-articular injection. Therefore the mixing or simultaneous administration with other intra-articular injectables is not recommended.

Storage

Store below 25 °C.

Do not freeze.

Protect from light ☼.

Sterile product for single use only ☑.

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