


ANTIBACTERIAL DRESSING

MEDIHONEY®

Antibacterial Honey Tulle Dressing
Antibakterieller Honig-Tüllverband
Pansement en tulle au miel antibactérien
Miel Antibacteriana Apósito de Tul
Gazeforbinding med Antibakteriel Honning
Antibakteerinen lääkehunajavoideside
Antibakteriellt honungsomslag av tyll
Medicazione di tulle al miele antibatterico
Αντιβακτηριακό Μέλι Επίδεσμος Γάζας
Antibacteriële Honing Wondverband
Compressa de Gaze com Mel Bactericida
Antibakteriální medový obvaz z tylu
گاز پانسومان عسل ضد باکتری
Antibakteriyel Ballı Tül Pansuman
ضمادة التول بالعسل المضاد للبكتيريا
Antibakterijska obveza iz tkanine z medom
Gasbandasje med antibakteriell honning
Opatrunek tiulowy z miodem antybakteryjnym

Antibacterial Honey Gel Sheet
Antibakterieller Honig-Gelverband
Gel Pansement au Miel Antibactérien
Miel Antibacteriana Lámina de Gel
Gelforbinding med Antibakteriel Honning
Antibakteerinen hunajageelityyny
Antibakteriell honung Gelförband
Foglio gel al miele antibatterico
Αντιβακτηριακό Μέλι Φύλλο Γέλης
Antibacteriële Honing Gelkompres
Lámina de Gel com Mel Bactericida
Antibakteriální medový Gel Sheet
ژل ورقه ای عسل ضد باکتری
Antibakteriyel Ballı Yaprak Jel
الهلام بالعسل المضاد للبكتيريا
ضمادة شريحة
Antibakterijska obloga z gelom
Gelkompres med antibakteriell honning
Arkus z żelowy z miodem antybakteryjnym

Antibacterial Honey Apinate™ Dressing
Antibakterieller Honig-Alginatverband
Pansement Alginate au Miel Antibactérien
Miel Antibacteriana Apósito de Alginato
Alginatbandage med Antibakteriel Honning
Antibakteerinen hunaja-alginaattisidos
Antibakteriell honung Alginatförband
Medicazione con alginato al miele antibatterico
Αντιβακτηριακό Μέλι Αλγινικό Επίθεμα
Antibacteriële Honing Alginaat Wondverband
Penso de Alginato com Mel Bactericida
Antibakteriální medový alginátový obvaz
آپینات™ یرتکاب دض لسع نامسناپ
Antibakteriyel Ballı Alginat Pansuman
آپینات™ بالعسل المضاد للبكتيريا
ضمادة ألجینات
Antibakterijska obloga z alginatom
Antibakteriell alginatbandasje med honning
Opatrunek Apinate™ z miodem antybakteryjnym

Antibacterial Honey Apinate™ Dressing product components from the UK and New Zealand. Packaged and Sterilized in Canada.

Antibacterial Honey Tulle Dressing product components from New Zealand. Packaged and Sterilized in Canada

Antibacterial Honey Gel Sheet made in the USA.

 Derma Sciences, Inc.
 104 Shorting Road
 Toronto, Ontario M1S 3S4
 USA Int: +1 609 514 4744
 Canada Int: +1 416 299 4003

Derma Sciences Europe Ltd.
 P.O. Box 220
 5 High Street
 Maidenhead SL6 1JN UK
 Int: +44 1628 779 461
 Freephone (UK only): 0800 071 3912

www.dermasciences.com

EC REP

Scanlan Group B.V.
 Postbus 75664,
 1118 ZS, Schiphol Triport
 The Netherlands,
 Int: +31 20 653-0553

Australian Sponsor
 Emergo Australia
 201 Sussex Street
 Darling Park, Tower II, Level 20
 Sydney, NSW 2000 Australia

CE
 0120

STERILE | R

It's not just any honey – it's MediHoney®

All of the products in the Derma Sciences® MediHoney® range contain MediHoney® Antibacterial Honey, which is a standardized medical grade of antibacterial honey collected from *Leptospermum* sp. (Manuka). This honey has been processed using advanced scientific knowledge to maximize the product's potential to heal.

Antibacterial Honey Tulle Dressing

MediHoney® Antibacterial Honey Tulle Dressing (Tulle Dressing) is a sterile, non-adherent, non-absorbent wound dressing comprising of a strong woven dressing impregnated with MediHoney® Antibacterial Honey.

Antibacterial Honey Apinate™ Dressing

MediHoney® Antibacterial Honey Apinate™ Dressing (Apinate™ Dressing) is a sterile, non-adherent wound dressing comprising of calcium alginate fibre impregnated with MediHoney® Antibacterial Honey. The dressing will absorb wound exudate to form a soft gel.

Antibacterial Honey Gel Sheet

MediHoney® Antibacterial Honey Gel Sheet (Gel Sheet Dressing) is a sterile, non-adherent wound dressing comprising of MediHoney® Antibacterial Honey and sodium alginate. These ingredients are combined using patented technology to create a honey dressing with a unique texture. The dressing will form a soft gel when in contact with wound exudate. Not recommended for dry wounds.

Actions:

- The MediHoney® Tulle, Apinate™ and Gel Sheet Dressings are for use on acute and chronic wounds. Their actions include:
 - assisting wound healing (providing a moist wound healing environment)
 - protecting the wound by creating an antibacterial barrier against bacteria including antibiotic resistant strains
 - fast, effective autolytic debridement to clean the wound, removing slough and debris
 - rapidly removing malodour
 - reducing the risk of infection from bacteria

Indications:

	Antibacterial Honey Tulle Dressing	Antibacterial Honey Apinate™ Dressing	Antibacterial Honey Gel Sheet
Leg / foot ulcers	✓	✓	✓
Pressure ulcers	✓	✓	✓
Infected wounds	✓	✓	✓
Sloughy wounds	✓	✓	✓
Malodorous wounds	✓	✓	✓
Diabetic wounds	✓	✓	✓
Donor and recipient graft sites	✓	✓	✓
Burns	✓	✓	✓
Superficial wounds such as cuts, scratches, abrasions	✓	✓	✓
Surgical wounds	✓	✓	✓
Necrotic wounds	✓	✓	✓

Serious wounds should be managed under the supervision of a health care professional.

Directions for use:

1. If required protect the wound edges with MediHoney® Barrier Cream.
 2. The Gel Sheet Dressing has protective liners on either side of the dressing. Remove these prior to use.
 3. Place the dressing onto the wound surface ensuring that the dressing is in contact with the wound bed. Dressings can be placed alongside or on top of one another to cover large or deep areas. Alternatively the wound may be filled with MediHoney® Antibacterial Medical Honey® before applying the Dressing. The Dressing may be cut with sterile scissors or folded to the shape of the wound.
 4. A suitable sterile secondary dressing should be applied to cover the MediHoney® Dressing. The secondary dressing should be sufficiently absorbent to manage wound exudate.
 5. The wound may be cleaned between dressing changes in accordance with standard wound cleaning protocol.
- Suitable for use under compression dressings.

Frequency of change:

- The MediHoney® Dressings require changing when they have been diluted by wound exudate and/or when the dressing has gelled. The dressings may be left in place for up to 7 days depending on the amount of wound exudate.
- All MediHoney® Dressings are non-adherent and therefore easy to remove without disrupting the wound bed and new tissue.
- Heavily exuding wounds require close monitoring and may require daily dressing changes.

Contraindications:

Do not use on individuals who have a known sensitivity to any component present in the dressing.

Observations:

The wound may appear larger in size and depth as dead tissue is cleaned away. This increase should be followed by an improved appearance of the wound.

Precautions:

- Due to the properties of honey some patients may experience some discomfort such as a stinging or drawing sensation upon application. This will subside over time. If discomfort is of an unacceptable level remove the dressing and wash the affected area.
- When used on dry wounds it may be helpful to moisten the dressing to aid removal.
- Use with caution on heavily bleeding wounds
- Do not use if pack is damaged.
- Product is only sterile until pack is opened.

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