

X-HA[®] volume[®]

HYALURONIC ACID VOLUMIZING IMPLANT TO REDEFINE FACIAL VOLUMES

INDICATIONS

X-HA[®] Volume is indicated to restore volumes of:

- Cheek
- Cheekbones
- Chin
- Face oval
- Temporal area

FORMULATION

Its viscosity is twice as high as X-HA^{3®}

- Resorbable monophasic implant
- 23 mg/ml highly cross-linked hyaluronic acid
- High molecular weigh (2,5 millions Da)

INJECTION

X-HA[®] Volume can be injected in deep dermis

ADVANTAGES

/ HIGH VISCOSITY

X-HA Volume[®] has a very high viscosity to resist to deformations and create volume.

X-HA Volume[®] is twice as volumizing as X-HA^{3®}


/ PERFECT HOMOGENEITY

Easy injection

Thanks to its high viscosity, X-HA[®] Volume can be easily injected with a non traumatic needle.

/ EASY USE

X-HA Volume[®] is packed in a box of 2 1ml-syrines to better adapt volume to the necessary corrections: 1 to 3ml if needed.



Injected quantities*

| | 1ml | 1,5ml | 2ml |
|------------------|------|-------|-----|
| Cheekbones | 25% | 8% | 67% |
| Cheek | 40% | | 60% |
| Temporal area | 100% | | |
| Nasolabial folds | 33% | | 67% |
| Ovale | 80% | | 20% |

* Observational multicenter study on 47 patients