



Information for use

LUITPOLD PHARMA GMBH · Zielstattstraße 9 · D-81379 München

Mobilat[®] Gel

for the percutaneous treatment of
contusions, bruises, muscular strains,
effusions, painful movements.

Composition

adrenocortical extract	1.0 g
mucopolysaccharide polysulfuric acid ester	0.2 g
salicylic acid	2.0 g
gel base	to 100.0 g

The active ingredients of Mobilat Gel rapidly and completely penetrate the skin, inhibit inflammation, and reduce effusions.

Mobilat Gel increases local blood and lymph flow, activates the metabolism and the regenerative forces of the afflicted connective tissue region.

Pain and swellings diminish rapidly after application of Mobilat Gel. It is well tolerated and has a pleasantly cooling effect.

Directions for use

Unless otherwise directed by the physician Mobilat Gel is applied several times daily. It should not be brought into contact with open wounds or mucous membranes.

Presentation

Tubes of 50 g of gel



LUITPOLD
SANKYO GROUP

engl./arab.

0067 10 000 001 0-07/95